

Antibodies Product Datasheet

Product Name:	Anti-COVID-19 N protein mAb IgM (S06FC)
Catalog Number:	CA0005
Reactivity:	COVID-19
Applications:	ELISA ECL Immunogold
Dilutions:	ELISA 1:5000-1:100000, WB 1:1000-1:10000
Modification:	Not Modification Specific
Source:	Humanized mAb IgM
Purification:	The antibody was affinity-purified from 293T expression by affinity-chromatography using specific immunogen.
Immunogen:	The recombinant SARS-CoV-2 (2019-nCoV) N Protein
Concentration:	
Gene Name:	SARS-CoV-2 N, 2019-nCoV N, COVID-19 N
Protein Name:	SARS-CoV-2 Nucleoprotein , 2019-nCoV Nucleoprotein, COVID-19 Nucleoprotein
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS
Specificity:	This antibody detects SARS-CoV-2 (2019-nCoV) Nucleoprotein
Storage and Stability:	-20°C/1 year