

Antibodies Product Datasheet

Product Name:	E2F-1 (Acetyl-Lys120) Antibody
Catalog Number:	D12058
Reactivity:	Human:K120, Mouse:K115, Rat:K118
Applications:	WB, ELISA
Dilutions:	WB: 1:500-10000 ELISA: 1:10000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized acetyl-peptide from human protein at AA range: 100-170
Concentration:	1 mg/ml
Gene Name:	E2F1 RBBP3
Protein Name:	E2F transcription factor 1
Other Name:	E2F1 RBBP3
Human Swiss Prot No.:	Q01094
Formulation:	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous E2F-1 (Acetyl-Lys120) protein
Storage and Stability:	-20°C/1 year