

Antibodies Product Datasheet

Product Name:	E2F-1 (Acetyl-K125) Antibody
Catalog Number:	D12059
Reactivity:	Human:K125, Mouse:K120, Rat:K123
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized Acetyl peptide derived from human E2F-1. at AA range: K125
Concentration:	1 mg/ml
Gene Name:	E2F1 RBBP3
Protein Name:	E2F-1
Other Name:	Transcription factor E2F1 (E2F-1) (PBR3) (Retinoblastoma-associated protein 1) (RBAP-1) (Retinoblastoma-binding protein 3) (RBBP-3) (pRB-binding protein E2F-1)
Human Swiss Prot No.:	Q01094
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of E2F-1 (Acetyl-K125), It doesn't react with total protein.
Storage and Stability:	-20°C/1 year