

Antibodies Product Datasheet

Product Name:	E2F-1 Antibody
Catalog Number:	D12059N
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 WB 1:500-2000, ELISA 1:10000-20000 IHC 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 100-170
Concentration:	1 mg/ml
Gene Name:	E2F1 RBBP3
Protein Name:	E2F transcription factor 1
Other Name:	E2F1 RBBP3
Human Swiss Prot No.:	Q01094
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous E2F-1 protein
Storage and Stability:	-20°C/1 year