

Antibodies Product Datasheet

Product Name:	E2F-4 (Acetyl-Lys96) Antibody
Catalog Number:	D12062
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized Acetyl-peptide derived from human E2F4 around the Acetylation site of Lys96. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	E2F4
Protein Name:	Transcription factor E2F4
Other Name:	E2F4; Transcription factor E2F4; E2F-4
Human Swiss Prot No.:	Q16254
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Acetyl-E2F-4 (K96) Polyclonal Antibody detects endogenous levels of E2F-4 protein only when acetylated at K96.
Storage and Stability:	-20°C/1 year