

Antibodies Product Datasheet

Product Name:	CRM1 (Acetyl-Lys568) Antibody
Catalog Number:	D12069
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized acetyl-peptide derived from the Internal region of human CRM1 around the acetylation site of K568.
Concentration:	1 mg/ml
Gene Name:	XPO1
Protein Name:	Exportin-1
Other Name:	XPO1; CRM1; Exportin-1; Exp1; Chromosome region maintenance 1 protein homolog
Human Swiss Prot No.:	O14980
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Acetyl-CRM1 (K568) Polyclonal Antibody detects endogenous levels of CRM1 protein only when acetylation at K568.
Storage and Stability:	-20°C/1 year