

Antibodies Product Datasheet

Product Name:	HNF-3 α /3 β /3 γ (Acetyl-Lys264/253/211) Antibody
Catalog Number:	D12070
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized Acetyl-peptide derived from human FOXA1 around the Acetylation site of Lys265. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	FOXA1/FOXA2/FOXA3
Protein Name:	Hepatocyte nuclear factor 3-alpha/Hepatocyte nuclear factor 3-beta/Hepatocyte nuclear factor 3-gamma
Other Name:	FOXA1; HNF3A; TCF3A; Hepatocyte nuclear factor 3-alpha; HNF-3-alpha; HNF-3A; Forkhead box protein A1; Transcription factor 3A; TCF-3A; FOXA2; HNF3B; TCF3B; Hepatocyte nuclear factor 3-beta; HNF-3-beta; HNF-3B; Forkhead box protein A2; Transcription factor
Human Swiss Prot No.:	P55317
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Acetyl-HNF-3 α /3 β /3 γ (K264/253/211) Polyclonal Antibody detects endogenous levels of HNF-3 α /3 β /3 γ protein only when acetylated at K264/253/211.
Storage and Stability:	-20°C/1 year