

Antibodies Product Datasheet

Product Name:	GATA-2/3 (Acetyl-Lys336/304) Antibody
Catalog Number:	D12078
Reactivity:	Human:K336/304, Mouse:K336/304, Rat:K336/305
Applications:	WB, ELISA
Dilutions:	WB: 1:500-10000 ELISA: 1:10000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized acetyl-peptide from human protein at AA range: 320-390
Concentration:	1 mg/ml
Gene Name:	GATA2/GATA3
Protein Name:	GATA binding protein 2/GATA binding protein 3
Other Name:	GATA2
Human Swiss Prot No.:	P23769/P23771
Formulation:	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous GATA-2/3 (Acetyl-Lys336/304) protein
Storage and Stability:	-20°C/1 year