

Antibodies Product Datasheet

Product Name:	Histone H2A.Z Antibody
Catalog Number:	D12083N
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 1-30
Concentration:	1 mg/ml
Gene Name:	H2AFZ H2AZ
Protein Name:	Histone H2A.Z (H2A/z)
Other Name:	Histone H2A.Z (H2A/z)
Human Swiss Prot No.:	P0C0S5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous Histone H2A.Z
Storage and Stability:	-20°C/1 year