

Antibodies Product Datasheet

Product Name:	Histone H2B (Acetyl-Lys24/25) Antibody
Catalog Number:	D12091
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Acetyl peptide from human protein at AA range: 24/25
Concentration:	1 mg/ml
Gene Name:	Histone H2B
Protein Name:	Histone H2B
Other Name:	H2BFS; Histone H2B type F-S; Histone H2B.s; H2B/s;H2BK2425AC
Human Swiss Prot No.:	Q96A08/P33778/P62807
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous Histone H2B when Acetyl occurs at Lys24 or 25)
Storage and Stability:	-20°C/1 year