

Antibodies Product Datasheet

Product Name:	HMG-1 Antibody
Catalog Number:	D12124N
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human HMG-1.
Concentration:	1 mg/ml
Gene Name:	HMGB1
Protein Name:	High mobility group protein B1
Other Name:	HMGB1; HMG1; High mobility group protein B1; High mobility group protein 1; HMG-1
Human Swiss Prot No.:	P09429
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HMG-1 Polyclonal Antibody detects endogenous levels of HMG-1 protein.
Storage and Stability:	-20°C/1 year