

Antibodies Product Datasheet

Product Name:	Karyopherin α 2 Antibody
Catalog Number:	D12140N
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300. Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human Karyopherin α 2.
Concentration:	1 mg/ml
Gene Name:	KPNA2
Protein Name:	Importin subunit alpha-2
Other Name:	KPNA2; RCH1; SRP1; Importin subunit alpha-2; Karyopherin subunit alpha-2; RAG cohort protein 1; SRP1-alpha
Human Swiss Prot No.:	P52292
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Karyopherin α 2 Polyclonal Antibody detects endogenous levels of Karyopherin α 2 protein.
Storage and Stability:	-20°C/1 year