

Antibodies Product Datasheet

Product Name:	NFkB-p65 (Acetyl-Lys218) Antibody
Catalog Number:	D12160
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized Acetyl-peptide derived from human NFkB-p65 around the Acetylation site of Lys218. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	RELA
Protein Name:	Transcription factor p65
Other Name:	RELA; NFKB3; Transcription factor p65; Nuclear factor NF-kappa-B p65 subunit; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3
Human Swiss Prot No.:	Q04206
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Acetyl-NFkB-p65 (K218) Polyclonal Antibody detects endogenous levels of NFkB-p65 protein only when acetylated at K218.
Storage and Stability:	-20°C/1 year