

Antibodies Product Datasheet

Product Name:	EF-1 α 1/2 (Acetyl-Lys146) Antibody
Catalog Number:	D12217
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EF-1 α 1/2 (Acetyl Lys146)
Concentration:	1 mg/ml
Gene Name:	EEF1A1 EEF1A EF1A LENG7
Protein Name:	EF-1 α 1/2 (Acetyl Lys146)
Other Name:	Elongation factor 1-alpha 1 (EF-1-alpha-1) (Elongation factor Tu) (EF-Tu) (Eukaryotic elongation factor 1 A-1) (eEF1A-1) (Leukocyte receptor cluster member 7)
Human Swiss Prot No.:	P68104/Q05639/Q5VTE0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat EF-1 α 1/2 (Acetyl Lys146)
Storage and Stability:	-20°C/1 year