

Antibodies Product Datasheet

Product Name:	Huntingtin (Acetyl-Lys442) Antibody
Catalog Number:	D12221
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Huntingtin (Acetyl Lys442)
Concentration:	1 mg/ml
Gene Name:	HTT HD IT15
Protein Name:	Huntingtin (Acetyl Lys442)
Other Name:	Huntingtin (Huntington disease protein) (HD protein)
Human Swiss Prot No.:	P42858
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Huntingtin (Acetyl Lys442)
Storage and Stability:	-20°C/1 year