

Antibodies Product Datasheet

Product Name:	ACSS2 (Acetyl-Lys418) Antibody
Catalog Number:	D12223
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ACSS2 (Acetyl Lys418)
Concentration:	1 mg/ml
Gene Name:	ACSS2 ACAS2
Protein Name:	ACSS2 (Acetyl Lys418)
Other Name:	Acetyl-coenzyme A synthetase, cytoplasmic (EC 6.2.1.1) (Acetate--CoA ligase) (Acetyl-CoA synthetase) (ACS) (AceCS) (Acyl-CoA synthetase short-chain family member 2) (Acyl-activating enzyme)
Human Swiss Prot No.:	Q9NR19
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse ACSS2 (Acetyl Lys418)
Storage and Stability:	-20°C/1 year