

Antibodies Product Datasheet

Product Name:	Histone H1 (Acetyl-Lys66/63/64) Antibody
Catalog Number:	D12225
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Histone H1 (Acetyl Lys66/63/64)
Concentration:	1 mg/ml
Gene Name:	HIST1H1A H1F1
Protein Name:	Histone H1 (Acetyl Lys66/63/64)
Other Name:	Histone H1.1 (Histone H1a)
Human Swiss Prot No.:	Q02539/P16403/P16402/P10412
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Histone H1 (Acetyl Lys66/63/64)
Storage and Stability:	-20°C/1 year