

Antibodies Product Datasheet

Product Name:	Histone H2B (Acetyl-Lys35) Antibody
Catalog Number:	D12226
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Histone H2B (Acetyl Lys35)
Concentration:	1 mg/ml
Gene Name:	HIST1H2BC H2BFL; HIST1H2BE H2BFH; HIST1H2BF H2BFG; HIST1H2BG H2BFA; HIST1H2BI H2BFK
Protein Name:	Histone H2B (Acetyl Lys35)
Other Name:	Histone H2B type 1-C/E/F/G/I (Histone H2B.1 A) (Histone H2B.a) (H2B/a) (Histone H2B.g) (H2B/g) (Histone H2B.h) (H2B/h) (Histone H2B.k) (H2B/k) (Histone H2B.l) (H2B/l)
Human Swiss Prot No.:	P62807/P58876/Q93079/O60814/Q99880/Q99879/Q99877/Q5QNW6/P57053
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Histone H2B (Acetyl Lys35)
Storage and Stability:	-20°C/1 year