

Antibodies Product Datasheet

Product Name:	HMG-I/HMG-Y (Acetyl-Lys71) Antibody
Catalog Number:	D12229
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HMG-I/HMG-Y (Acetyl Lys71)
Concentration:	1 mg/ml
Gene Name:	HMGA1 HMG1Y
Protein Name:	HMG-I/HMG-Y (Acetyl Lys71)
Other Name:	High mobility group protein HMG-I/HMG-Y (HMG-I(Y)) (High mobility group AT-hook protein 1) (High mobility group protein A1) (High mobility group protein R)
Human Swiss Prot No.:	P17096
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat HMG-I/HMG-Y (Acetyl Lys71)
Storage and Stability:	-20°C/1 year