

Antibodies Product Datasheet

Product Name:	MYH-pan (Acetyl-Lys1394) Antibody
Catalog Number:	D12230
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MYH-pan (Acetyl Lys1394)
Concentration:	1 mg/ml
Gene Name:	MYH1
Protein Name:	MYH-pan (Acetyl Lys1394)
Other Name:	Myosin-1 (Myosin heavy chain 1) (Myosin heavy chain 2x) (MyHC-2x) (Myosin heavy chain IIx/d) (MyHC-IIx/d) (Myosin heavy chain, skeletal muscle, adult 1)
Human Swiss Prot No.:	P12882/Q9UKX2/P11055/Q9Y623/P13533/P12883/P13535
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat MYH-pan (Acetyl Lys1394)
Storage and Stability:	-20°C/1 year