

Antibodies Product Datasheet

Product Name:	GR (Acetyl-Lys494) Antibody
Catalog Number:	D12240
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GR (Acetyl Lys494)
Concentration:	1 mg/ml
Gene Name:	NR3C1 GRL
Protein Name:	GR (Acetyl Lys494)
Other Name:	Glucocorticoid receptor (GR) (Nuclear receptor subfamily 3 group C member 1)
Human Swiss Prot No.:	P04150
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat GR (Acetyl Lys494)
Storage and Stability:	-20°C/1 year