

Antibodies Product Datasheet

Product Name:	NCoA-3 (Acetyl-Lys616) Antibody
Catalog Number:	D12241
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NCoA-3 (Acetyl Lys616)
Concentration:	1 mg/ml
Gene Name:	NCOA3 AIB1 BHLHE42 RAC3 TRAM1
Protein Name:	NCoA-3 (Acetyl Lys616)
Other Name:	Nuclear receptor coactivator 3 (NCoA-3) (EC 2.3.1.48) (ACTR) (Amplified in breast cancer 1 protein) (AIB-1) (CBP-interacting protein) (pCIP) (Class E basic helix-loop-helix protein 42) (bHLHe42) (Receptor-associated coactivator 3) (RAC-3) (Steroid recepto
Human Swiss Prot No.:	Q9Y6Q9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human NCoA-3 (Acetyl Lys616)
Storage and Stability:	-20°C/1 year