

Antibodies Product Datasheet

Product Name:	Actin-pan (Acetyl-Lys63/61/62) Antibody
Catalog Number:	D12244
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Actin-pan (Acetyl Lys63/61/62)
Concentration:	1 mg/ml
Gene Name:	ACTA2 ACTSA ACTVS GIG46
Protein Name:	Actin-pan (Acetyl Lys63/61/62)
Other Name:	Actin, aortic smooth muscle (Alpha-actin-2) (Cell growth-inhibiting gene 46 protein)
Human Swiss Prot No.:	P62736/P60709/P68032/P63261/P63267/P68133
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Actin-pan (Acetyl Lys63/61/62)
Storage and Stability:	-20°C/1 year