

Antibodies Product Datasheet

Product Name:	Keratin-pan (Acetyl-Lys185) Antibody
Catalog Number:	D12245
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Keratin-pan (Acetyl Lys185)
Concentration:	1 mg/ml
Gene Name:	KRT2 KRT2A KRT2E
Protein Name:	Keratin-pan (Acetyl Lys185)
Other Name:	Keratin, type II cytoskeletal 2 epidermal (Cytokeratin-2e) (CK-2e) (Epithelial keratin-2e) (Keratin-2 epidermis) (Keratin-2e) (K2e) (Type-II keratin Kb2)
Human Swiss Prot No.:	P35908/Q01546/P12035/P13647/P02538/P04259/P48668/Q3SY84/Q14CN4/Q86Y46/Q7RTS7/O95678/Q5XKE5/P08729/P05787/Q9NSB2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Keratin-pan (Acetyl Lys185)
Storage and Stability:	-20°C/1 year