

Antibodies Product Datasheet

Product Name:	C/EBP β (Acetyl-Lys265) Antibody
Catalog Number:	D12253
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human C/EBP β (Acetyl Lys265)
Concentration:	1 mg/ml
Gene Name:	CEBPB LAP TCF5 PP9092
Protein Name:	C/EBP β (Acetyl Lys265)
Other Name:	CCAAT/enhancer-binding protein beta (C/EBP beta) (Liver activator protein) (Nuclear factor NF-IL6) (Transcription factor 5) (TCF-5)
Human Swiss Prot No.:	P17676
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat C/EBP β (Acetyl Lys265)
Storage and Stability:	-20°C/1 year