

Antibodies Product Datasheet

Product Name:	E2F-2 (Acetyl-Lys122) Antibody
Catalog Number:	D12258
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human E2F-2 (Acetyl Lys122)
Concentration:	1 mg/ml
Gene Name:	E2F2
Protein Name:	E2F-2 (Acetyl Lys122)
Other Name:	Transcription factor E2F2 (E2F-2)
Human Swiss Prot No.:	Q14209
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse E2F-2 (Acetyl Lys122)
Storage and Stability:	-20°C/1 year