

Antibodies Product Datasheet

Product Name:	EKLF (Acetyl-Lys274) Antibody
Catalog Number:	D12259
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EKLF (Acetyl Lys274)
Concentration:	1 mg/ml
Gene Name:	KLF1 EKLF
Protein Name:	EKLF (Acetyl Lys274)
Other Name:	Krueppel-like factor 1 (Erythroid krueppel-like transcription factor) (EKLF)
Human Swiss Prot No.:	Q13351
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human EKLF (Acetyl Lys274)
Storage and Stability:	-20°C/1 year