

Antibodies Product Datasheet

Product Name:	HIF-1 α (Acetyl-Lys532) Antibody
Catalog Number:	D12267
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HIF-1 α (Acetyl Lys532)
Concentration:	1 mg/ml
Gene Name:	HIF1A BHLHE78 MOP1 PASD8
Protein Name:	HIF-1 α (Acetyl Lys532)
Other Name:	Hypoxia-inducible factor 1-alpha (HIF-1-alpha) (HIF1-alpha) (ARNT-interacting protein) (Basic-helix-loop-helix-PAS protein MOP1) (Class E basic helix-loop-helix protein 78) (bHLHe78) (Member of PAS protein 1) (PAS domain-containing protein 8)
Human Swiss Prot No.:	Q16665
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat HIF-1 α (Acetyl Lys532)
Storage and Stability:	-20°C/1 year