

Antibodies Product Datasheet

Product Name:	Histone H2A (Acetyl-Lys119) Antibody
Catalog Number:	D12268
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Histone H2A (Acetyl Lys119)
Concentration:	1 mg/ml
Gene Name:	HIST1H2AG H2AFP; HIST1H2AI H2AFC; HIST1H2AK H2AFD; HIST1H2AL H2AFI; HIST1H2AM H2AFN
Protein Name:	Histone H2A (Acetyl Lys119)
Other Name:	Histone H2A type 1 (H2A.1) (Histone H2A/p)
Human Swiss Prot No.:	P0C0S8/P16104
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Histone H2A (Acetyl Lys119)
Storage and Stability:	-20°C/1 year