

Antibodies Product Datasheet

Product Name:	Histone H2A (Acetyl-Lys15) Antibody
Catalog Number:	D12270
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:500-1000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Histone H2A (Acetyl Lys15)
Concentration:	
Gene Name:	HIST1H2AG/HIST1H2AI/HIST1H2AK/HIST1H2AL/HIST1H2AM/HIST2H2AA3/HIST2H2AA4/HIST3H2A
Protein Name:	Histone H2A type 1/Histone H2A type 2/Histone H2A type 3
Other Name:	HIST1H2AG; H2AFP; HIST1H2AI; H2AFC; HIST1H2AK; H2AFD; HIST1H2AL; H2AFI; HIST1H2AM; H2AFN; Histone H2A type 1; H2A.1; Histone H2A/p; HIST2H2AA3; H2AFO; HIST2H2AA; HIST2H2AA4; Histone H2A type 2-A; Histone H2A.2; Histone H2A/o; HIST3H2A; Histone H2A type 3;
Human Swiss Prot No.:	P0C0S8/Q6FI13/Q7L7L0
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Histone H2A (Acetyl Lys15) protein.
Storage and Stability:	-20°C/1 year