

Antibodies Product Datasheet

Product Name:	Histone H2A.Z (Acetyl-Lys7) Antibody
Catalog Number:	D12279
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:500-1000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Histone H2A.Z (Acetyl Lys7)
Concentration:	
Gene Name:	H2AFZ
Protein Name:	Histone H2A.Z
Other Name:	H2AFZ; H2AZ; Histone H2A.Z; H2A/z; YM3299
Human Swiss Prot No.:	P0C0S5
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Histone H2A.Z (Acetyl Lys7) protein.
Storage and Stability:	-20°C/1 year