

Antibodies Product Datasheet

Product Name:	HMG-1 (Acetyl-Lys3) Antibody
Catalog Number:	D12308
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HMG-1 (Acetyl Lys3)
Concentration:	1 mg/ml
Gene Name:	HMGB1 HMG1
Protein Name:	HMG-1 (Acetyl Lys3)
Other Name:	High mobility group protein B1 (High mobility group protein 1) (HMG-1)
Human Swiss Prot No.:	P09429
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat HMG-1 (Acetyl Lys3)
Storage and Stability:	-20°C/1 year