

Antibodies Product Datasheet

Product Name:	KLF13 (Acetyl-Lys166) Antibody
Catalog Number:	D12313
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KLF13 (Acetyl Lys166)
Concentration:	1 mg/ml
Gene Name:	KLF13 BTEB3 NSLP1
Protein Name:	KLF13 (Acetyl Lys166)
Other Name:	Krueppel-like factor 13 (Basic transcription element-binding protein 3) (BTE-binding protein 3) (Novel Sp1-like zinc finger transcription factor 1) (RANTES factor of late activated T-lymphocytes 1) (RFLAT-1) (Transcription factor BTEB3) (Transcription fac
Human Swiss Prot No.:	Q9Y2Y9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse KLF13 (Acetyl Lys166)
Storage and Stability:	-20°C/1 year