

Antibodies Product Datasheet

Product Name:	Karyopherin α 2 (Acetyl-Lys22) Antibody
Catalog Number:	D12323
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Karyopherin α 2 (Acetyl Lys22)
Concentration:	1 mg/ml
Gene Name:	KPNA2 RCH1 SRP1
Protein Name:	Karyopherin α 2 (Acetyl Lys22)
Other Name:	Importin subunit alpha-2 (Karyopherin subunit alpha-2) (RAG cohort protein 1) (SRP1-alpha)
Human Swiss Prot No.:	P52292
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Karyopherin α 2 (Acetyl Lys22)
Storage and Stability:	-20°C/1 year