

Antibodies Product Datasheet

Product Name:	Histone H1 (Di Methyl-Lys25) Antibody
Catalog Number:	D12839
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:500-1000
Modification:	Meythl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Histone H1 (Di Methyl Lys25)
Concentration:	
Gene Name:	HIST1H1A/HIST1H1C/HIST1H1D/HIST1H1E
Protein Name:	Histone H1.1/Histone H1.2/Histone H1.3/Histone H1.4
Other Name:	H1K25ME2; HIST1H1A; H1F1; Histone H1.1; Histone H1a; HIST1H1C; H1F2; Histone H1.2; Histone H1c; Histone H1d; Histone H1s-1; HIST1H1D; H1F3; Histone H1.3; Histone H1c; Histone H1s-2; HIST1H1E; H1F4; Histone H1.4; Histone H1b; Histone H1s-4; YM3286
Human Swiss Prot No.:	Q02539/P16403/P16402/P10412
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Histone H1 (Di Methyl Lys25) protein.
Storage and Stability:	-20°C/1 year