

Antibodies Product Datasheet

Product Name:	Histone H2B (Tri Methyl-Lys5) Antibody
Catalog Number:	D12844
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:1000-2000
Modification:	Meythl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Histone H2B (Tri Methyl Lys5)
Concentration:	
Gene Name:	HIST1H2BC
Protein Name:	Histone H2B type 1-A/Histone H2B type 1-B/Histone H2B type 1-C/E/F/G/I
Other Name:	H2BK5ME3; HIST1H2BA; TSH2B; Histone H2B type 1-A; Histone H2B, testis; Testis-specific histone H2B; HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B.f; H2B/f; HIST1H2BC; H2BFL; HIST1H2BE; H2BFH; HIST1H2BF; H2BFG; HIST1H2BG; H2BFA; HIST1H
Human Swiss Prot No.:	Q96A08/P33778/P62807
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Histone H2B (Tri Methyl Lys5) protein.
Storage and Stability:	-20°C/1 year