

Antibodies Product Datasheet

Product Name:	C3aR Antibody
Catalog Number:	G028
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human C3AR1. AA range:175-224
Concentration:	1 mg/ml
Gene Name:	C3AR1
Protein Name:	C3a anaphylatoxin chemotactic receptor
Other Name:	C3AR1; AZ3B; C3R1; HNFAG09; C3a anaphylatoxin chemotactic receptor; C3AR; C3a-R
Human Swiss Prot No.:	Q16581
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	C3aR Polyclonal Antibody detects endogenous levels of C3aR protein.
Storage and Stability:	-20°C/1 year