

Antibodies Product Datasheet

Product Name:	AR- β 3 Antibody
Catalog Number:	G033
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ADRB3. AA range:250-299
Concentration:	1 mg/ml
Gene Name:	ADRB3
Protein Name:	Beta-3 adrenergic receptor
Other Name:	ADRB3; ADRB3R; B3AR; Beta-3 adrenergic receptor; Beta-3 adrenoreceptor; Beta-3 adrenoceptor
Human Swiss Prot No.:	P13945
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AR- β 3 Polyclonal Antibody detects endogenous levels of AR- β 3 protein.
Storage and Stability:	-20°C/1 year