

Antibodies Product Datasheet

Product Name:	CELSR3 Antibody
Catalog Number:	G058
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CELSR3. AA range:3195-3244
Concentration:	1 mg/ml
Gene Name:	CELSR3
Protein Name:	Cadherin EGF LAG seven-pass G-type receptor 3
Other Name:	CELSR3; CDHF11; EGFL1; FMI1; KIAA0812; MEGF2; Cadherin EGF LAG seven-pass G-type receptor 3; Cadherin family member 11; Epidermal growth factor-like protein 1; EGF-like protein 1; Flamingo homolog 1; hFmi1; Multiple epidermal growth factor-Q9NYQ7
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CELSR3 Polyclonal Antibody detects endogenous levels of CELSR3 protein.
Storage and Stability:	-20°C/1 year