

Antibodies Product Datasheet

Product Name:	CXCR-7 Antibody
Catalog Number:	G065
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CXCR7. AA range:311-360
Concentration:	1 mg/ml
Gene Name:	CXCR7
Protein Name:	C-X-C chemokine receptor type 7
Other Name:	CXCR7; CMKOR1; GPR159; RDC1; C-X-C chemokine receptor type 7; CXC-R7; CXCR-7; Chemokine orphan receptor 1; G-protein coupled receptor 159; G-protein coupled receptor RDC1 homolog; RDC-1
Human Swiss Prot No.:	P25106
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CXCR-7 Polyclonal Antibody detects endogenous levels of CXCR-7 protein.
Storage and Stability:	-20°C/1 year