

Antibodies Product Datasheet

Product Name:	CRF-RI Antibody
Catalog Number:	G072
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CRHR1. AA range:75-124
Concentration:	1 mg/ml
Gene Name:	CRHR1
Protein Name:	Corticotropin-releasing factor receptor 1
Other Name:	CRHR1; CRFR; CRFR1; CRHR; Corticotropin-releasing factor receptor 1; CRF-R-1; CRF-R1; CRFR-1; Corticotropin-releasing hormone receptor 1; CRH-R-1; CRH-R1
Human Swiss Prot No.:	P34998
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CRF-RI Polyclonal Antibody detects endogenous levels of CRF-RI protein.
Storage and Stability:	-20°C/1 year