

Antibodies Product Datasheet

Product Name:	EDG-2 Antibody
Catalog Number:	G083
Reactivity:	Human, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EDG2. AA range:5-54
Concentration:	1 mg/ml
Gene Name:	LPAR1
Protein Name:	Lysophosphatidic acid receptor 1
Other Name:	LPAR1; EDG2; LPA1; Lysophosphatidic acid receptor 1; LPA receptor 1; LPA-1; Lysophosphatidic acid receptor Edg-2
Human Swiss Prot No.:	Q92633
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EDG-2 Polyclonal Antibody detects endogenous levels of EDG-2 protein.
Storage and Stability:	-20°C/1 year