

Antibodies Product Datasheet

Product Name:	EMR1 Antibody
Catalog Number:	G090
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EMR1. AA range:818-867
Concentration:	1 mg/ml
Gene Name:	EMR1
Protein Name:	EGF-like module-containing mucin-like hormone receptor-like 1
Other Name:	EMR1; TM7LN3; EGF-like module-containing mucin-like hormone receptor-like 1; EGF-like module receptor 1; EMR1 hormone receptor
Human Swiss Prot No.:	Q14246
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EMR1 Polyclonal Antibody detects endogenous levels of EMR1 protein.
Storage and Stability:	-20°C/1 year