

Antibodies Product Datasheet

Product Name:	ETL Antibody
Catalog Number:	G100
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ELTD1. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	ELTD1
Protein Name:	EGF latrophilin and seven transmembrane domain-containing protein 1
Other Name:	ELTD1; ETL; EGF; latrophilin and seven transmembrane domain-containing protein 1; EGF-TM7-latrophilin-related protein; ETL protein Q9HBW9
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ETL Polyclonal Antibody detects endogenous levels of ETL protein.
Storage and Stability:	-20°C/1 year