

Antibodies Product Datasheet

Product Name:	FSHR Antibody
Catalog Number:	G101
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IF 1:100-300 ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FSHR. AA range:278-327
Concentration:	1 mg/ml
Gene Name:	FSHR
Protein Name:	Follicle-stimulating hormone receptor
Other Name:	FSHR; LGR1; Follicle-stimulating hormone receptor; FSH-R; Follitropin receptor
Human Swiss Prot No.:	P23945
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FSHR Polyclonal Antibody detects endogenous levels of FSHR protein.
Storage and Stability:	-20°C/1 year