

Antibodies Product Datasheet

Product Name:	FPR1 Antibody
Catalog Number:	G102
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FPR1. AA range:155-204
Concentration:	1 mg/ml
Gene Name:	FPR1
Protein Name:	fMet-Leu-Phe receptor
Other Name:	FPR1; fMet-Leu-Phe receptor; fMLP receptor; N-formyl peptide receptor; FPR; N-formylpeptide chemoattractant receptor
Human Swiss Prot No.:	P21462
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FPR1 Polyclonal Antibody detects endogenous levels of FPR1 protein.
Storage and Stability:	-20°C/1 year