

Antibodies Product Datasheet

Product Name:	Frizzled-9 Antibody
Catalog Number:	G114
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FZD9. AA range:542-591
Concentration:	1 mg/ml
Gene Name:	FZD9
Protein Name:	Frizzled-9
Other Name:	FZD9; FZD3; Frizzled-9; Fz-9; hFz9; FzE6; CD antigen CD349
Human Swiss Prot No.:	O00144
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Frizzled-9 Polyclonal Antibody detects endogenous levels of Frizzled-9 protein.
Storage and Stability:	-20°C/1 year