

Antibodies Product Datasheet

Product Name:	G2A Antibody
Catalog Number:	G130
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR132. AA range:293-342
Concentration:	1 mg/ml
Gene Name:	GPR132
Protein Name:	Probable G-protein coupled receptor 132
Other Name:	GPR132; G2A; Probable G-protein coupled receptor 132; G2 accumulation protein
Human Swiss Prot No.:	Q9UNW8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	G2A Polyclonal Antibody detects endogenous levels of G2A protein.
Storage and Stability:	-20°C/1 year