

Antibodies Product Datasheet

Product Name:	GPR175 Antibody
Catalog Number:	G146
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR175. AA range:198-247
Concentration:	1 mg/ml
Gene Name:	TPRA1
Protein Name:	Transmembrane protein adipocyte-associated 1
Other Name:	TPRA1; GPR175; PP6566; Transmembrane protein adipocyte-associated 1; Integral membrane protein GPR175
Human Swiss Prot No.:	Q86W33
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR175 Polyclonal Antibody detects endogenous levels of GPR175 protein.
Storage and Stability:	-20°C/1 year