

Antibodies Product Datasheet

Product Name:	GPR18 Antibody
Catalog Number:	G147
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR18. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	GPR18
Protein Name:	N-arachidonyl glycine receptor
Other Name:	GPR18; GPCRW; N-arachidonyl glycine receptor; NAGly receptor; G-protein coupled receptor 18
Human Swiss Prot No.:	Q14330
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR18 Polyclonal Antibody detects endogenous levels of GPR18 protein.
Storage and Stability:	-20°C/1 year